SARS-CoV-2 Antigen Detection Kit (Colloidal Gold Method)

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Nr.	Name	Probennamen
1	Influenza A H1N1	Influenza A H1N1 (New Cal/20/99) Culture Fluid (1 mL)
2	Influenza A H1N1 (Seasonal)	Influenza A H1N1 (Seasonal) Culture Fluid (1 mL)
3	Influenza A H3N2	Influenza A H3N2 (Brisbane/10/07) Culture Fluid (1 mL)
4	Influenza A H5N1	Influenza A H5N1 Culture Fluid (1 mL)
5	Influenza B Yamagata	Influenza B (Yamagata/16/88) Culture Fluid (1 mL)
6	Influenza B Victoria	Influenza B (Victoria/2/87) Culture Fluid (1 mL)
7	Haemophilus Influenzae	Haemophilus Influenzae Culture Fluid (1 mL)
8	Respiratory Syncytial Virus Type A	Respiratory Syncytial Virus Type A (RSV-A) Culture Fluid (1 mL)
9	Respiratory Syncytial Virus Type B	Respiratory Syncytial Virus Type B (RSV-B) Culture Fluid (1 mL)
10	Rhinovirus Type 1A	Rhinovirus Type 1A Culture Fluid (1 mL)
11	Rhinovirus B14	Rhinovirus B14 Culture Fluid (1 mL)
12	Adenovirus Type 1	Adenovirus Type 1 Culture Fluid (1 mL)
13	Adenovirus Type 2	Adenovirus Type 2 Culture Fluid (1 mL)
14	Adenovirus Type 3	Adenovirus Type 3 Culture Fluid (1 mL)
15	Adenovirus Type 4	Adenovirus Type 4 Culture Fluid (1 mL)
16	Adenovirus Type 5	Adenovirus Type 5 Culture Fluid (1 mL)
17	Adenovirus Type 7A	Adenovirus Type 7A Culture Fluid (1 mL)
18	Epstein-Barr Virus	Epstein-Barr Virus (EBV) Culture Fluid (1 mL)
19	Measles Virus	Measles Virus Culture Fluid (1 mL)
20	Cytomegalovirus	Cytomegalovirus (CMV) (Strain: AD-169) Culture Fluid (1 mL)

Nr.	Name	Probennamen
21	Rotavirus	Rotavirus Culture Fluid (1 mL)
22	Norovirus	Norovirus Group I (Recombinant) Culture Fluid (1 mL)
23	Enterovirus A71	Enterovirus A71 Culture Fluid (1 mL)
24	Mycobacterium Tuberculosis	Mycobacterium Tuberculosis Culture Fluid (1 mL)
25	Mumps Virus	Mumps Virus (Isolate 1) Culture Fluid (1 mL)
26	Streptococcus Pyogenes	Streptococcus Pyogenes Culture Fluid (1 mL)
27	Streptococcus Pneumoniae	Streptococcus Pneumoniae Culture Fluid (1 mL)
28	Candida Albicans	Candida Albicans Culture Fluid (1 mL)
29	Varicella Zoster Virus	Varicella Zoster Virus (VZV) Strain: 82 Culture Fluid (1 mL)
30	Mycoplasma Pneumoniae	Mycoplasma Pneumoniae M129, Titered (1 MI)
31	Chlamydia Pneumoniae	Chlamydia Pneumoniae, Titered (1 MI)
32	Bordetella Pertussis	Bordetella Pertussis Culture Fluid (1 MI)
33	Legionella Pneumophila	Legionella Pneumophila Culture Fluid (1 MI)
34	Staphylococcus Aureus	Staphylococcus Aureus Culture Fluid (1 MI)
35	Staphylococcus Epidermidis	Staphylococcus Epidermidis Culture Fluid (1 MI)
36	Candida Albicans	Candida Albicans Culture Fluid (1 MI)
37	Parainfluenza Virus 1	Parainfluenza Virus 1 Culture Fluid (1 MI)
38	Parainfluenza Virus 2	Parainfluenza Virus 2 Culture Fluid (1 MI)
39	Parainfluenza Virus 3	Parainfluenza Virus 3 Culture Fluid (1 MI)
40	Parainfluenza Virus 4	Parainfluenza Virus 4 Culture Fluid (1 MI)
41	Coronavirus HKU1	Coronavirus HKU1 Culture Fluid (1 mL)
42	Coronavirus OC43	Coronavirus OC43 Culture Fluid (1 mL)
43	Coronavirus NL63	Coronavirus NL63 Culture Fluid (1 mL)
44	Coronavirus 229E	Coronavirus 229E Culture Fluid (1 mL)
45	Coronavirus MERS-Cov	MERS-Cov Culture Fluid (1 mL)

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Literaturverzeichnis

- World Health Organization: Clinical management of severe acute respiratory infection when SARS-CoV-2 (nCoV) infection is suspected: Interim Guidance. 12 January, 2020.
- Chen, Yu, Lanjuan Li, et al. „SARS-CoV-2: virus dynamics and host response.“ The Lancet Infectious Diseases, 2020, Vol. 20(5): 515-516.

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Ningbo Lvtang Biotechnology Co., Ltd.
No.9 Dongpu Road, Chengdong Industrial Area,
Daxu Town, Xiangshan County, Ningbo City,
Zhejiang Province, China
Isa.yuan@lvtangbio.com
+86-574-65002553

EC REP Lotus N L B.V.
Koningin Julianaplein 10, 1e Verd
2595AA, The Hague, Niederlande
peter@lotusnl.com

Nasentupfer

Ningbo HLS Medical Products Co., Ltd.
No.3918-3928, Zhouxi Road, Zhouxiang Town,
Cixi City, Zhejiang 315324, China

EC REP MedPath GmbH
Mies-van-der-Rohe-Strasse 8
80807 München, Deutschland
Email: info@medpath.pro

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SARS-CoV-2 Antigen Detection Kit (Colloidal Gold Method)

For Self-testing in vitro diagnostic use only

Instruction for Use

Packaging Specification

1 Test / Bag, 1 Test / Box, 2 Tests / Box, 5 Tests / Box, 20 Tests / Box, 25 Tests / Box
Model: CovAg-N

Intended Use

This product is used for qualitative detection of SARS-CoV-2 antigen in human nasal swab samples. For those with a history of exposure, but no symptoms and symptoms of coronavirus pneumonia disease, screening should be performed at the initial stage of symptoms (0 to 7 days) after infection. The SARS-CoV-2 virus belongs to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. The test is dedicated for symptomatic individuals up to 7 days after symptoms occurrence, or asymptomatic people with known exposure or after suspected contact might be tested. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. This product provides preliminary test results. A negative result does not preclude infection with a SARS-CoV-2 virus and is not used as the sole basis for treatment or other management decisions. The kit is intended for layperson use and for people age 5 years or above. Children should be assisted by an adult. People over 65 years of age should seek help with the test.

Test Principle

This kit uses immunochromatography for detection. Colloidal gold was used as a marker to qualitatively detect SARS-CoV-2 N protein antigen in human nasal swab samples. The specimen will move forward along the test card under capillary action. If the specimen contains a SARS-CoV-2 antigen, the antigen will bind to the colloidal gold-labeled SARS-CoV-2 monoclonal antibody to form an immune complex. The immune complex will be fixed on the membrane, and SARS-CoV-2 monoclonal antibody will be captured, form the fuchsia line, display will be SARS-CoV-2 antigen positive; If the line does not show color, the negative result will be displayed. The test card also contains a quality control line C, which shall appear fuchsia color regardless of whether there is a detection line.

Main Components

1 Test / Bag, 1 Test / Box: 1 piece of all components
20 Tests/Box: 20 pieces of all components
2 Tests/Box: 2 pieces of all components
25 Tests/Box: 25 pieces of all components
5 Tests/Box: 5 pieces of all components

Precaution

- For self-testing in vitro diagnostic use only. Do not use after expiration date.
- Read the instructions carefully before using the kit, and strictly control the reaction time. If you do not follow the instructions, you will get inaccurate results.
- Guard against moisture, do not open the aluminum platinum bag before it is ready for testing. Do not use the aluminum foil bag when it is damaged or the test card is damp. Do not take the aluminium foil bag desiccant orally.
- Balance all reagents and specimens to room temperature (10°C-30°C) before use.
- Do not replace the components in this kit with components in other kits.
- Do not dilute the specimen for testing, otherwise you may get inaccurate results.
- The kit shall be stored in strict accordance with the conditions specified in this manual. Please do not store the kit under freezing conditions.
- The test methods and results must be interpreted in strict accordance with this specification.
- Negative results will occur with this kit if the SARS-CoV-2 antigen titer in the specimen falls below the minimum.

Storage Conditions & Period of Validity


- Store at 2°C-30°C, and it is valid for 24 months.
- After the aluminum foil bag is unsealed, the test card should be used as soon as possible within 1 hour.

Preparatory Steps

- Disinfect the surface where you will open the test kit. Remove and lay out all contents of the test kit on a clean, flat surface.
- Please blow your nose and clear the nostril before taking the test.
- Wash hands with soap and water for at least 20 seconds. If soap and water are not available, use hand sanitizer. Dry your hands completely before taking the test.
- You need to prepare a timer or any device with timer function.

Implementation

Please read the manual carefully before using it. The samples to be tested, detection reagents and other materials for detection shall be balanced to room temperature (10-30°C), and the test shall be conducted at room temperature (10-30°C). Please note, that the test card must be used within 1 hour once the foil pouch has been opened.

- Collect the Nasal Swab Sample. Use the sterile swab provided in the kit, carefully insert the swab into one nostril of the user. The swab tip should be inserted up to 1.5-2 cm until resistance is met.
- Roll the swab 5 times along the mucosa inside the nostril to ensure that both mucus and cells are collected. Use the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities.
- Withdraw the swab from the nasal cavity.
- Peel the aluminum foil, slightly tilt the extraction buffer tube. Stir the swab at least 5 times.
- Move the swab up and down 3 times while squeezing the extraction buffer tube in order to get as many specimens as possible from the swab. Take out the swab and put on the cap.
- Open the aluminium foil bag and take out the test card. Apply 2-3 drops of extracted specimen to the specimen hole (S) of the test card.
- Read and interpret the test result after 15-20 minutes. The test result should not be read and interpreted after 30 minutes.
- Put all the components into the specimen bag and disposed with of household waste to prevent pollution of the environment and wash your hands for at least 20s.

Interpretation of Test Results

Negative Result: if there is only a quality control line C and the detection line T is colorless, indicating that SARS-CoV-2 antigen has not been detected and the result is negative.

A negative test result does not rule out possibility of infection, an infection may be present even if your test result is negative. You should continue to comply with all applicable rules regarding contact with others and protective measures. In case of suspicion, please repeat the test after 1-2 days, as the novel coronavirus may not be detected accurately in all phases of an infection.

Positive Result: if both the quality control line C and the detection line T appear, SARS-CoV-2 antigen has been detected and the result is positive for antigen.

Look very closely! The T line can be very faint. A positive test result suggests there is currently a suspicion of COVID-19 infection, if your test result is positive, you should immediately contact a doctor/family doctor or the local health authority and ask for their professional opinions, you should comply with local self-isolation guidelines and have a COVID-19 nucleic acid PCR confirmatory test carried out to confirm the infection.

Invalid Result: if the quality control line C is not observed, it will be invalid regardless of whether there is a detection line T (as shown in the figure on the left), and the test shall be conducted again.

An invalid test result is possibly caused by faulty test execution, please repeat the test; If test result is still invalid, contact a doctor or COVID-19 test center for their professional opinions and immediately contact the manufacturer or local supplier.

Quality Control

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Limitations of Inspection Methods

- This kit is for the detection of human nasal swab samples extracts only. Results from other specimens may be erroneous.
- The accuracy of the test depends on the process of collecting the sample, and improper sample collection, improper sample storage, or repeated freezing and thawing of the sample can affect the test results.
- This product can only qualitatively detect SARS-CoV-2 N protein antigen in human nasal swab samples and cannot accurately determine the content of antigen in samples.
- The detection results of this reagent are only for clinical reference and shall not be used as the only basis for clinical diagnosis and treatment. The clinical management of patients shall be comprehensively considered in combination with their symptoms/signs, medical history, other laboratory tests and treatment reactions.
- Due to the limitation of antigen detection reagent methodology, for negative detection results, it is recommended to use nucleic acid detection. Please seek follow-up care with your physician or healthcare provider as additional testing may be necessary.
- For individuals without COVID-19 symptoms and/or individuals who live in areas with low numbers of COVID-19 infections and without known exposure to COVID-19, more false positive results may occur. Testing of individuals without symptoms should be limited to contacts of confirmed or probable cases or to other epidemiological reasons to suspect a COVID-19 infection and should be followed by additional confirmatory testing with a molecular test. Please seek follow-up care with your physician or healthcare provider.
- This reagent is a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time, based on the doctor's

diagnosis.

8. Analysis of the possibility of false negative results:

- 1) Unreasonable sample collection, transportation and treatment.
- 2) There was no virus in the sample, or the virus content in the sample was very low, which was lower than the critical detection concentration of the reagent, and the low-concentration sample could not be detected. Viral genetic variants may lead to changes in antigenic determinants, resulting in false negative results.
- 3) The optimal sample type and the optimal sampling time after infection (peak virus) are not verified, therefore, it is possible to avoid false negatives by collecting samples at multiple sites in multiple fractions from the same person.
- 4) The virus in the sample may have been inactivated. Although nucleic acid fragments may remain, viral antigens may have been destroyed and inactivated.

Performance

1. Minimum detection limit

The LOD and repeatability of SARS-CoV-2 strain or recombinant SARS-CoV-2 nucleocapsid protein samples studied, and the positive detection rate of SARS-CoV-2 Antigen Detection Kit(Colloidal Gold Method) is $\geq 90\%$, the repeatability tests all meet the requirements.

According to the test results in Table 1 Test results of SARS-CoV-2 Culture Fluid (Heat Inactivated) the limit of detection(LOD) of SARS-CoV-2 Antigen Detection Kit (Colloidal Gold Method), the limit of detection(LOD) concentration is $1.17 \times 10^3 \text{TCID}_{50}/\text{mL}$

According to the test results in Table 2 Test results of recombinant SARS-CoV-2 nucleocapsid protein the limit of detection(LOD) of SARS-CoV-2 Antigen Detection Kit(Colloidal Gold Method), the limit of detection(LOD) concentration is $1.0 \times 10^3 \text{pg}/\text{mL}$.

2. Sensitivity, Specificity and Accuracy Nasal Swab

A total of 446 various types of samples are tested. Compared with the virellaSARS-CoV-2 seqc real time RT-PCR Kit, the SARS-CoV-2 Antigen Detection Kit (Colloidal Gold Method) See below Table for details shows Sensitivity:95.50% (95% CI: 90.40%-97.90%), Specificity: 98.70% (95% CI: 96.80%-99.50%), Accuracy: 97.80% (95% CI:95.90%-98.80%).

Clinical samples	virellaSARS-CoV-2 seqc real time RT-PCR Kit		Total
	Positive	Negative	
SARS-CoV-2 Antigen Detection Kit (Colloidal Gold Method)	Positive 126	Negative 4	130
	Negative 6	310	316
Total	132	314	446

Sensitivity: 95.50% (95% CI: 90.40%-97.90%)

Specificity: 98.70% (95% CI: 96.80%-99.50%)

Accuracy: 97.80% (95% CI: 95.90%-98.80%)

3. Hook effect

Testing of $3.5 \times 10^3 \text{TCID}_{50}/\text{mL}$ SARS-CoV-2 virus after heat inactivation resulted in no Hook effect for the assay.

4. Cross-reactivity

Cross-reactivity with the following organisms was investigated. Samples that were positive for the following organisms tested negative with this assay. After heat inactivation of $3.5 \times 10^3 \text{TCID}_{50}/\text{mL}$ SARS-CoV-2 virus, the following microorganisms (Table below) showed no cross reaction with the detection reagent at the following concentrations. Source: zeptomatrix

No	Name	Sample name
1	Influenza A H1N1	Influenza A H1N1 (New Cal/20/99) Culture Fluid (1 mL)
2	Influenza A H1N1 (Seasonal)	Influenza A H1N1 (Seasonal) Culture Fluid (1 mL)
3	Influenza A H3N2	Influenza A H3N2 (Brisbane/10/07) Culture Fluid (1 mL)
4	Influenza A H5N1	Influenza A H5N1 Culture Fluid (1 mL)
5	Influenza B Yamagata	Influenza B (Yamagata/16/89) Culture Fluid (1 mL)
6	Influenza B Victoria	Influenza B (Victoria/20/07) Culture Fluid (1 mL)
7	Haemophilus influenzae	Haemophilus influenzae Culture Fluid (1 mL)
8	Respiratory Syncytial Virus Type A	Respiratory Syncytial Virus Type A (RSV-A) Culture Fluid (1 mL)
9	Respiratory Syncytial Virus Type B	Respiratory Syncytial Virus Type B (RSV-B) Culture Fluid (1 mL)
10	Rhinovirus Type 1A	Rhinovirus Type 1A Culture Fluid (1 mL)
11	Rhinovirus B14	Rhinovirus B14 Culture Fluid (1 mL)
12	Adenovirus Type 1	Adenovirus Type 1 Culture Fluid (1 mL)
13	Adenovirus Type 2	Adenovirus Type 2 Culture Fluid (1 mL)
14	Adenovirus Type 3	Adenovirus Type 3 Culture Fluid (1 mL)
15	Adenovirus Type 4	Adenovirus Type 4 Culture Fluid (1 mL)
16	Adenovirus Type 5	Adenovirus Type 5 Culture Fluid (1 mL)
17	Adenovirus Type 7A	Adenovirus Type 7A Culture Fluid (1 mL)
18	Epstein-Barr Virus	Epstein-Barr Virus (EBV) Culture Fluid (1 mL)
19	Measles Virus	Measles Virus Culture Fluid (1 mL)
20	Cytomegalovirus	Cytomegalovirus (CMV) (Strain: AD-169) Culture Fluid (1 mL)
21	Rotavirus	Rotavirus Culture Fluid (1 mL)
22	Norovirus	Norovirus Group 1 (Recombinant) Culture Fluid (1 mL)
23	Enterovirus A71	Enterovirus A71 Culture Fluid (1 mL)
24	Mycobacterium Tuberculosis	Mycobacterium Tuberculosis Culture Fluid (1 mL)
25	Mumps Virus	Mumps Virus (Isolate 1) Culture Fluid (1 mL)
26	Streptococcus Pyogenes	Streptococcus Pyogenes Culture Fluid (1 mL)
27	Streptococcus Pneumoniae	Streptococcus Pneumoniae Culture Fluid (1 mL)
28	Candida Albicans	Candida Albicans Culture Fluid (1 mL)
29	Varicella Zoster Virus	Varicella Zoster Virus (VZV) Strain: 82 Culture Fluid (1 mL)
30	Mycoplasma Pneumoniae	Mycoplasma Pneumoniae M129, Titered (1 MI)
31	Chlamydia Pneumoniae	Chlamydia Pneumoniae, Titered (1 MI)
32	Bordetella Pertussis	Bordetella Pertussis Culture Fluid (1 MI)
33	Legionella Pneumophila	Legionella Pneumophila Culture Fluid (1 MI)
34	Staphylococcus Aureus	Staphylococcus Aureus Culture Fluid (1 MI)
35	Staphylococcus Epidermidis	Staphylococcus Epidermidis Culture Fluid (1 MI)
36	Candida Albicans	Candida Albicans Culture Fluid (1 MI)
37	Parainfluenza Virus 1	Parainfluenza Virus 1 Culture Fluid (1 MI)
38	Parainfluenza Virus 2	Parainfluenza Virus 2 Culture Fluid (1 MI)
39	Parainfluenza Virus 3	Parainfluenza Virus 3 Culture Fluid (1 MI)
40	Parainfluenza Virus 4	Parainfluenza Virus 4 Culture Fluid (1 MI)
41	Coronavirus HKU1	Coronavirus HKU1 Culture Fluid (1 mL)
42	Coronavirus OC43	Coronavirus OC43 Culture Fluid (1 mL)
43	Coronavirus NL63	Coronavirus NL63 Culture Fluid (1 mL)
44	Coronavirus 229E	Coronavirus 229E Culture Fluid (1 mL)
45	Coronavirus MERS-Cov	MERS-Cov Culture Fluid (1 mL)

5. Interfering Substances

The following substances (Table below), which are naturally present in respiratory specimens or may be artificially introduced into the nasal cavity, nasopharynx, were evaluated at the following concentrations when tested with the assay and found not to affect the performance of the test.

Drug classification	Drug name	Dosage	
Allergic symptom relief medication	Histamine Dihydrochloride	0.1mg/time	
	Antiviral drugs	Interferon a	0.5ml/time
		Zanamivir	10mg/time
		Ribavirin	0.3g/time
		Oseltamivir	75mg/time
		Parainfluenza	300mg/time
		Lopinavir	200mg/time
		Ritonavir	50mg/time
		Arbidol	0.2g/time
		Antibiotics	Levofloxacin
Azithromycin			1g/time
Ceftriaxone	1g/time		
Meropenem	200mg/time		
Systemic antibacterial drug	Tobramycin	0.1ml/time	
	Other drugs	Phenylephrine	1g/time
		Oxymetazoline	1g/time
		0.9% Sodium Chloride	1g/time
		Natural Soothing Alkali	1g/time
		Beclomethasone	1g/time
		Dexamethasone	1g/time
		Haloperidol	1g/time
		Triamcinolone	1g/time
		Budesonide	1g/time
		Mometasone	1g/time
		Fluticasone	1g/time
		Fluticasone Propionate	1g/time

Bibliography

1. World Health Organization: Clinical management of severe acute respiratory infection when SARS-CoV-2 (nCoV) infection is suspected: Interim Guidance. 12 January, 2020.
2. Chen, Yu, Lanjuan Li, et al. „SARS-CoV-2: virus dynamics and host response.“ The Lancet Infectious Diseases, 2020, Vol. 20(5): 515-516.

Index of Symbols

In vitro diagnostic medical device	Do not reuse	Use-by date
Caution: Please refer to the instruction in the packaging	Caution: Please refer to the instruction in the packaging	Manufacturer
Temperature limit	Batch number	Authorized representative in the European Community
Keep dry.	Keep away from sun light	Do not use if package is damaged
Date of manufacture	Contains sufficient for <n> tests	CE mark

Ningbo Lvtang Biotechnology Co., Ltd.
No.9 Dongqiu Road, Chengdong Industrial Area,
Daxu Town, Xiangshan County, Ningbo City,
Zhejiang Province, China
lsa.yuan@lvtangbio.com
+86-574-85002553

EC REP Lotus N L.B.V.
Koningin Julianaplein 10, 1e Verd
2595AA, The Hague, Netherlands.
peter@lotusn.com

Disposable Swabs

Ningbo HLS Medical Products Co., Ltd.
No.3818-3928, Zhouxi Road, Zhouxiang Town,
Cixi City, Zhejiang 315324, China

EC REP MedPath GmbH
Mies-van-der-Rohe-Strasse 8
80807 Munich, Germany
Email:info@medpath.pro

Approval Date and Revision Date of the Instruction for use

Approval Date: May, 15th, 2020

Revision Date: Feb. 1st, 2022

Date of Issue: Feb. 5th, 2022