

# SARS-CoV-2 **Antigen Rapid Test**

(Colloidal Gold Immunochromatography)

For Self-testing





Nantong Egens Biotechnology Co., Ltd

Email: Egens@Egens-bio.cn Web address: www.egens-bio.cn

2022 Edition I

## **SARS-CoV-2** Antigen Rapid Test

### for Self-testing



#### **Product Features**

- Specimen type:Nasal Aspirate Fluid/Nasal Swab/
- Easy to collect samples, simple operation, without professional equipment and laboratory
- Testing time:15-20 minutes
- convenient transportation
- Applicable to self- testing at home

#### **Sensitivity And Specificity**

	RT-PCR		
SARS-CoV-2 Antigen Rapid Test	Positive	Negative	TOTAL
Positive	205	1	206
Negative	14	599	613
Total	219	600	819

**Sensitivity**=205/219×100%=93.61% (Wilson 95%CI: 89.51% - 96.46%)

**Specificity**=599/600×100%=99.83% (Wilson 95%CI: 99.07% - 100.00%)

**Total accuracy**=(205+599)/819×100%=98.17%

(Wilson 95%CI: 97.00% - 98.97%)

#### **Kit Components**

1 Test/Box	5 Tests/Box	25 Tests/Box	
1 Test Cassette	5 Test Cassettes	25 Test Cassettes	
1 Sterilized swab	5 Sterilized swabs	25 Sterilized swabs	
1 Extraction Tube	5 Extraction	25 Extraction Tubes	
with Buffer	Tubes with Buffer	with Buffer	
1 Biosafety bag	5 Biosafety bags	25 Biosafety bags	
1 Instructions for	1 Instructions for	1 Instructions for	
use	use	use	
1 cardboard rack	1 cardboard rack	1 cardboard rack	











## SARS-CoV-2 Antigen Rapid Test for Self-testing


#### **Specimen Requirements and Test Procedure**



Wash your hands before starting







Check all components for possible damage and that all materials supplied are complete. Lay all the supplied materials on a clean, dry and flat surface.



Pinch the upper of the extraction tube and shake it more times to make buffer drop to the bottom of the tube completely. Tear off the aluminum foil seal of the extraction tube containing the buffer solution carefully. (\*Note: Keep your distance away with your face)







Insert the extraction tube into the circular hole of cardboard rack.



Take out the test cassette from the sealed foil pouch and place it on a clean and level surface.






Please remove excess mucus in the nasal cavity before sampling by swab.

Unpack the swab package, being careful NOT to touch the soft end, which is the absorbent tip.




Gently insert the swab into your nasal cavity 1-1.5cm deep (generally 1cm for children aged 2-14 years). Do not inserting deeper, avoid hurting your nostril. Slowly rotate the swab in the inside walls of nostril 5 times or more to collect more nasal secretions. Repeat steps in the other nostril using the same swab.



Insert the swab into the extraction tube with buffer, gently pick up the tube from cardboard rack, and then rotate the swab for 10~15 times to mix the collected specimen with the buffer as much as possible.





of the extraction tube, pinch the swab head with one hand across the outer wall of the extraction tube, and bend the swab to break it gently with another hand.



Close the cap of the extraction tube and press firmly.






Add 2 drops of the treated sample vertically into the sample well of the test cassette. And start the timer.



Read the result 15-20 minutes after adding the sample. Result got after 20 minutes is invalid.












All used test components should be put into the biosafety bag, sealed and discarded. (According to local garbage classification regulations)

\*Note: After completing all steps, wash hands with hand sanitizer.

### Interpretation of result







#### **Packing information**

Product Code	D0101TE	D0501EE	D2501AT
Format	1 Test/box,400 Box/Ctn	5 Test/box,200 Box/Ctn	25 Test/box,40 Box/Ctn
Box Size	13.5*7*2CM	15*7*4CM	19.5*13.2*8.5CM
Carton Size	57*38*43CM	63*38*43CM	55.5*41*45.5CM
G.W.	12.0KGS	16.0KGS	14.0KGS

### SARS-CoV-2 Antigen Rapid Test for Self-testing



EC Certificate No. 1434-IVDD-231/2022 EC Design-examination Directive 98/79/EC concerning in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies that manufactured by:

Nantong Egens Biotechnology Co., Ltd. Building 15, Building 12 (west), No. 1692 Xinghu Avenue, Nantong Economy & Technology Development Zone, 226010 Nantong, PEOPLE'S REPUBLIC OF CHINA.

> in vitro diagnostic medical devices for self-testing

#### SARS-CoV-2 Antigen Rapid Test Kit

The list of medical devices covered by this certificate is provided in the annex 1

in terms of design documentation, comply with requirements of Annex III (Section 6) to Directive 98/79/EC (as amended) implemented into Polish law, as evidenced by the audit conducted by the PCBC Validity of the Certificate: From 25.05.2022 to 27.05.2025

The date of issue of the Certificate: 25.05.2022 The date of the first issue of the Certificate: 25.05.2022



der the Contract No. MD-171/2021 on No: 332/2021 e bears the qualified signature.

Director
Medical Device Certification
Department

POLISH CENTRE FOR TESTING AND CERTIFICATION 02-844 Warsaw, 469 Pulawska Street, tel. +48 22 46 45 200, e-mail:pcbc@pcbc.gov.pl



#### ANNEX 1 TO THE CERTIFICATE VALID ONLY WITH CERTIFICATE

No. 1434-IVDD-231/2022

List of medical devices covered by the certificate: 1 Test/Box - REF. D0101TE 5 Tests/Box - REF. D0501EE 25 Tests/ Box - REF. D2501AT



Certificate bears the qualified signature. Warsaw, 24/05/2022

POLISH CENTRE FOR TESTING AND CERTIFICATION 02-844 Warsaw, 469 Pulawska Street, tel. +48 22 46 45 200, e-mail:pcbc@pcbc.gov.i

